



Contains Natural Rubber Latex

PowerFlex® AFD

ABSORBENT FOAM DRESSING



U.S. Patent No. 5,762,623
European Patent No. 0839222

PowerFlex® AFD is an all-in-one foam dressing and cohesive bandage that delivers instant application mid-game. A highly absorbent foam pad is attached to strong, flexible PowerFlex for quick and easy application. A waterproof film layer between the PowerFlex and foam pad prevents leakage and contamination. Fast, effective treatment for wounds, lacerations and turf burns. PowerFlex AFD consolidates three steps of typical treatment into one easy step!

Uses for PowerFlex® AFD:

- ✓ Cover lacerations
- ✓ Cover turf burns
- ✓ Cover abrasions


- Quick, all-in-one wound management system
- Excellent Blood Bandage - stops bleeding quickly
- Keeps you in the game!
- Made in the USA

Catalog #	Size	Rolls/Case	Colors
3722XX-016	2" x 2.5 yds.	16	○●
			○ (S) = sterile product/white only
3744XX-008	4" x 2.5 yds.	8	○●
			○ (S) = sterile product/white only

Available in White (WH) or Black (BK)

When ordering, please indicate catalog # and color code. For example, 3722WH-016 (XX = color code)

PowerFlex AFD Construction



- ➊ PowerFlex Layer - Strong & flexible, sticks to itself for quick & easy wrapping
- ➋ Waterproof Film - Prevents leakage and contamination
- ➌ Absorbent Foam - Foam pad draws in fluid while keeping wound site clean