



Contains
Natural
Rubber
Latex

PowerPro™

SYNTHETIC ADHESIVE TAPE



PowerPro™ is an adhesive athletic tape made without cotton, giving it less than 1% stretch. It delivers the performance of PowerTape and VictoryTape but in a versatile adhesive that sticks to anything, not just PowerFlex. PowerPro hand-tears with ease and has exceptional “lay in” to conform beautifully around ankles. PowerPro remains supportive and stable throughout games and practices, as well as in adverse playing conditions.

Uses for PowerPro™:

- ✓ Replace traditional cloth athletic tape
- ✓ Provide support & stability during games

- Less than 1% stretch - maintains integrity
- Stable support & range of motion control
- Crisp, easy hand-tear
- Sticks to anything - PowerFlex, pre-wrap, or skin
- Cotton free - won't absorb sweat or water
- Made in the USA

Catalog #	Size	Rolls/Case	Colors
ACP165XX-1"	1" x 15 yds.	48	○●
ACP165XX-1.5"	1.5" x 15 yds.	32	○●
ACP165XX-2"	2" x 15 yds.	24	○●

Available in White (WH) or Black (BK)

When ordering, please indicate catalog # and color code.

For example, ACP165WH-1" (XX = color code)