



Contains Natural Rubber Latex

VictoryTape™

ECONOMICAL + EASY TO HAND-TEAR



Patent Pending

VictoryTape™ is an economical cohesive athletic tape with easy hand-tear technology that sticks only to itself or to PowerFlex. It contains less than 1% stretch, eliminating the “give or slip” that exists when using traditional cotton tape. Apply over PowerFlex for an ultimate taping system of support, conformability and comfort. The system is sweat and water resistant, remaining supportive and stable throughout games and practices, even in adverse playing conditions.

Uses for VictoryTape™:

- ✓ Replace traditional cloth athletic tape
- ✓ Use in combination with PowerFlex for the Ultimate Taping System

- Easy hand-tear cohesive tape that sticks to itself or PowerFlex, not skin or hair
- Less than 1% stretch - eliminates the “give or slip” that is found with cloth tape
- Sweat and water resistant
- All-day performance
- Made in the USA

Catalog #	Size	Rolls/Case	Colors
ACP120XX-1"	1" x 15 yds.	48	○●
ACP120XX-1.5"	1.5" x 15 yds.	32	○●
ACP120XX-2"	2" x 15 yds.	24	○●

Available in White (WH) or Black (BK)

When ordering, please indicate catalog # and color code.
For example, ACP120WH-1" (XX = color code)

What's the difference between PowerTape® and VictoryTape™?

✓ PowerTape is the original and still gold standard of cohesive tapes with 41 lbs. of tensile strength and less than 1% stretch.

✓ VictoryTape is a more economical, easy to hand-tear version of PowerTape. It has 38 lbs. of tensile strength and less than 1% stretch.