






The Next Generation of Cohesive Bandages.



Quick-Stick Application in a Latex-Free Cohesive Bandage.

Features and Benefits:

- Quick-stick open cell foam clings to itself, stays in place and will not unravel
- Water-resistant cohesive is unaffected by water and holds up to wet dressings
- Soft and comfortable foam layer is non-irritating and gentle on skin
- Thin and lightweight for added comfort and prolonged wear
- Latex-free product

Applications:

A general purpose cohesive bandage for all applications including the treatment of edema, holding dressings in place, stabilizing IV's and as a light compression wrap.



Why is CoFlex LF² Better?

CoFlex LF² sticks better than other latex-free bandages.

CoFlex LF² is latex-free, but sticks like latex. The open cell foam on the surface of CoFlex LF² is designed to absorb and hold the latex-free cohesive. This allows for a quick-stick, long lasting cohesive bond.

CoFlex LF² is more comfortable than other cohesive bandages.

CoFlex LF² is designed with a thin, soft layer of specialty foam in the bandage. The foam provides a smooth surface that comes in contact with the skin, maintaining a comfort level that is unmatched by competitive products.

CoFlex LF² sticks over wet dressings.

CoFlex LF² was designed to stick over wet dressings. Traditionally latex-free bandages haven't performed well over wet dressings, forcing the industry to continue using latex bandages. CoFlex LF² is latex-free and its patented cellular foam technology enables it to stick and hold wet dressings in place.

CoFlex LF² is easier to apply and lasts longer than other cohesive bandages. A combination of Easy-Tear Patented Technology and Cellular Foam Patented Technology delivers a bandage that has easy hand-tear and will last until removed so there is no need to reapply the bandage due to unraveling.

U.S. Patent Nos. 5	ORDERING IN5,762,623 & 6,15	
Cellular Foam U.S	5. Patent Pending	
Catalog#	Size	Rolls/Case
9100XX-030	1" x 5 yds.	30
9150XX-048	1 ¹ /2" x 5 yds.	48
9200XX-036	2" x 5 yds.	36
9300XX-024	3" x 5 yds.	24
9400XX-018	4" x 5 yds.	18
9600XX-012	6" x 5 yds.	12
	Tan (TN), Red (RD), Green (NG), Purple (blorpack (CP)	
Sterile (S) = Ster	ile Product/Tan only	<i>י</i> .
9200(S)	2" x 5 yds.	24
9300(S)	3" x 5 yds.	24
9400(S)	4" x 5 yds.	18
9600(S)	6" x 5 yds.	12

CellularFoam EasyTear



Andover Healthcare, Inc.

9 Fanaras Drive, Salisbury, MA 01952

Toll Free: 800-432-6686 Phone: 978-465-0044 Fax: 978-462-0003 www.andoverhealthcare.com