

# PowerTape<sup>®</sup>

SELF-ADHERING SPORTS TAPE

LATEX  
Contains  
Natural  
Rubber  
Latex



Patent Pending

PowerTape<sup>®</sup> is a cohesive athletic tape that sticks only to itself or to PowerFlex. It contains less than 1% stretch, eliminating the “give or slip” that exists when using traditional cotton tape. With 41 lbs of tensile strength, PowerTape combined with PowerFlex create an ultimate taping system of maximum support, conformability and comfort. The system is sweat and water resistant, remaining supportive and stable throughout games and practices, even in adverse playing conditions.

## Uses for PowerTape<sup>®</sup>:

- ✓ Replace traditional cloth athletic tape
- ✓ Use in combination with PowerFlex for the Ultimate Taping System

- Hand-tear cohesive tape that sticks to itself or PowerFlex, not skin or hair
- Less than 1% stretch - eliminates the “give or slip” that is found with cloth tape
- Sweat and water resistant
- All-day performance
- Made in the USA

Catalog #	Size	Rolls/Case	Colors
ACP130XX-1"	1" x 15 yds.	48	○●
ACP130XX-1.5"	1.5" x 15 yds.	32	○●
ACP130XX-2"	2" x 15 yds.	24	○●

Available in White (WH) or Black (BK)

When ordering, please indicate catalog # and color code.

For example, ACP130WH-1" (XX = color code)

## PowerTape Application



### Step 1

Apply PowerFlex directly to the skin with proper unwind tension, removing all stretch.



### Step 2

Cover PowerFlex with PowerTape to provide the strongest taping system in existence.