

# CoFlex® AFD

ABSORBENT FOAM DRESSING



## All-in-One Foam Pad Dressing & Bandage

CoFlex® AFD's highly absorbent foam pad is attached to soft, quick-stick CoFlex® LF<sup>2</sup> for easy application. Between the LF<sup>2</sup> layer and the absorbent foam is a waterproof film layer that prevents leakage, drying out and contamination. Also available sterile.



### Key Features:

- ✓ Quick & easy to apply
- ✓ High absorption and retention of fluids (equivalent to 29 gauze pads)
- ✓ Flexible, controlled compression
- ✓ Made in the USA

### PRODUCT APPLICATIONS:

- Trauma Bandaging
- Post-op cranial surgery
- General wound dressing
- Cover/protect post surgery incisions
- Post bunion surgery

Catalog #	Size	Rolls/Case	Catalog #	Size	Rolls/Case	Colors
9720TN-016	2" x 2.5 yds.	16	9760TN-008	6" x 2.5 yds.	8	●
9720(S)						
9740TN-008	4" x 2.5 yds.	8	9840TN-006	4" x 7.5 yds.	6	●
9740(S)						

Choose From: Tan (TN), Sterile Tan (S)

When ordering, please indicate catalog # and color code  
 For example, 9720TN-016 (XX = color code)  
 Most products are available in bulk packaging, call for details.