

# Help prevent biting, tearing & chewing with:



## PowerFlex<sup>®</sup>

SELF-ADHERING ELASTIC BANDAGE

**NO CHEW**



- Strong - 23 lbs. tensile strength
- Straight hand tear
- Sweat and water resistant
- Controlled compression - will not constrict
- Quick, easy application - no scissors needed

Catalog #	Size	Rolls/Case
3840NC-018	4" x 5 yds.	18

Our bitter bandages are safe for you and your pet, however may leave a residue on your hands after use. Wearing gloves and washing hands thoroughly after application is recommended. NO CHEW PetFlex contains bitterants that have been tested on and found safe to use on domestic animals.

Andover Healthcare, Inc. | Phone: 800.432.6686 | [www.andoverhealthcare.com](http://www.andoverhealthcare.com)