

# CoFlex® AFD

ABSORBENT FOAM DRESSING



## All-in-One Foam Pad Dressing & Bandage

CoFlex® AFD's highly absorbent foam pad is attached to soft, quick-stick CoFlex® LF<sup>2</sup> for easy application. Between the LF<sup>2</sup> layer and the absorbent foam is a waterproof film layer that prevents leakage, drying out and contamination. Also available sterile.



### Key Features:

- ✓ Quick & easy to apply
- ✓ High absorption and retention of fluids (equivalent to 29 gauze pads)
- ✓ Non-adherent foam pad conforms but won't stick to wound bed
- ✓ Flexible, controlled compression
- ✓ Made in the USA

### PRODUCT APPLICATIONS:

- Trauma Bandaging
- Post-op cranial surgery
- General wound dressing
- Cover/protect post surgery incisions
- Post bunion surgery

| Catalog #  | Size          | Rolls/Case | Catalog #  | Size          | Rolls/Case | Colors |
|------------|---------------|------------|------------|---------------|------------|--------|
| 9720TN-016 | 2" x 2.5 yds. | 16         | 9760TN-008 | 6" x 2.5 yds. | 8          | ●      |
| 9720(S)    |               |            |            |               |            |        |
| 9740TN-008 | 4" x 2.5 yds. | 8          | 9840TN-006 | 4" x 7.5 yds. | 6          | ●      |
| 9740(S)    |               |            |            |               |            |        |

Choose From: Tan (TN), Sterile Tan (S)

When ordering, please indicate catalog # and color code  
 For example, 9720TN-016 (XX = color code)  
 Most products are available in bulk packaging, call for details.