

# CoFlex® + UMAFD

UNIVERSAL MEDIC ABSORBENT FOAM DRESSING



## All-in-One Trauma/Combat Bandage

CoFlex® UMAFD is a strong bandage with foam pad attached, plus an additional foam pad in the core for extra absorbency. It offers quick, easy application in high stress situations. It is a great alternative to existing combat bandages. Additionally, CoFlex® UMAFD has been heat tested at 130° > 90 days.

### Key Features:

- ✓ Quick & easy to apply
- ✓ High absorption and retention of fluids (equivalent to 29 gauze pads)
- ✓ Non-adherent foam pad conforms but won't stick to wound bed
- ✓ Flexible, controlled compression
- ✓ Able to withstand high temperatures
- ✓ Made in the USA

NO LATEX



Catalog #	Size	Rolls/Case	Colors
9340UM-018	4" x 3 yds.	18	●
9340UM-100	4" x 3 yds.	100	●

Available in Black only

Durable trauma bandage for situations requiring additional wound packing, strength for compression, and high heat/water resistance for any of the following:

- Military - Medics, Soldiers
- ER / Trauma Centers
- First Responders
- Commercial & Industrial Sites
- Emergency Medical Services
- First-Aid Kits

“

We were very happy with how easy the CoFlex AFD foam pad dressing was to apply. Patients we sent home with the product were able to change the dressings without difficulty and therefore more compliant. We used the CoFlex AFD foam pad dressing on many types of wounds and were pleased with the healing we saw with them all.

Amy Mahoney, Especialista em Tratamento de Feridas  
Boone Hospital, Wound Clinic  
Columbia, MO